


OPO

>> *Optical Parametric Oscillator [OPO] Package for deep tissue imaging* <<



- Excitation of red dyes
- Reduced photo toxicity and bleaching
- High penetration depth
- 1100 - 1600 nm @ 200 fs


Extended infrared multiphoton imaging

Most standard 2-photon microscopes are equipped with Ti:Sapphire lasers that deliver NIR laser radiation in the spectral range between 680-1040 nm. Therefore 2-photon microscopy is limited to blue, green and yellow dyes or fluorescent proteins - red dyes or red fluorescent proteins are almost excluded from this deep imaging technique. Exciting red dyes or red fluorescent proteins requires pulsed laser radiation in the far infrared (>1100 nm), which is delivered by OPOs.

Technique/physics behind the OPO

OPOs are passive laser-like devices that are designed for operation with mode-locked Ti:Sapphire lasers. The OPO converts the short Ti:Sa laser radiation into tuneable radiation of longer wavelengths. In 2-photon microscopy setups the Ti:Sapphire laser pumps the OPO at 740-880 nm (@ fixed wavelength) and the OPO delivers tunable fs laser radiation from 1100 to 1600 nm.

OPO setup



LaVision BioTec's flexible approach

LaVision BioTec's flexible approach LaVision BioTec's OPO package provides a flexible motorized beam splitter that allows to set the fraction of light used for pumping the OPO, while the other fraction excites blue/green/yellow fluorophores simultaneously. The fraction used for exciting blue/green/yellow fluorophores passes a negative pulse compressor to ensure laser pulse length < 150 fs within the sample. Ti:Sa laser and OPO wavelengths can be tuned independently to excite blue/green/yellow respectively red fluorophores most efficiently.

Benefits for 2-photon microscopy

Extended infrared radiation delivered by the OPO offers following advantages:

1. Access to red dyes/fluorescent protein
2. Significantly reduced photo toxicity and bleaching rate
3. Higher penetration depth (depending on the specimen up to a factor of 2)

Literature

1. 1419-24 *Biophysical J.* 2010 Feb; 98: 715–723. *Expanding Two-Photon Intravital Microscopy to the Infrared by Means of Optical Parametric Oscillator.* Herz J, Siffrin VD, Hauser AE, Brandt AU, Leuenberger T, Radbruch H, Zipp F, Niesner RA.

2. *Curr Opin Biotechnol.* 2009 Mar 25. *Infrared multiphoton microscopy: subcellular-resolved deep tissue imaging.* Andresen V, Alexander S, Heupel WM, Hirschberg M, Hoffman RM, Friedl P.

LaVision BioTec 

LaVision BioTec GmbH, Astastraße 14, D-33617 Bielefeld, Germany
Fon: +49 / 521 / 91 51 39-0, Fax: +49 / 521 / 91 51 39-10
Email: info@lavisvisionbiotec.com, Internet: www.lavisvisionbiotec.com